

FAQs V1.2

Last updated on: 20/12/2017
Inceptor Technologies Pvt. Ltd. (Innovision)
Firmware version: M01.12S01.04H09.00

1. Prospective User

- **What is the weight and size of Braille Me?**

Braille Me weighs just 580 g and its dimensions are 18.8cm x 13.8cm x 2.9cm.

- **What are the contents of the box when I first receive my Braille Me?**

The package contains the Braille Me device, a printed user guide, a braille user guide, Standard-A to Micro-B USB Cable, AC Charger and an 8GB SD card pre-inserted into the device.

- **How fast is the braille refresh rate?**

The refresh rate for the device stands at 300 milliseconds, which is just as quick as the time required by the human eye for a blink.

- **How much sound does it make while typing and refreshing braille?**

The sound generated while typing is pretty much similar to any computer keyboard. Refreshing the braille itself hardly makes a soft noise and can easily be used in a classroom environment.

- **What are the features of the stand-alone mode?**

In this mode, the device displays the content of files stored on an SD card. You can access previously prepared files from the SD card or create new files and save them for further use.

- **What is the maximum size of SD card supported? What file formats are supported?**

The device supports a card with a maximum storage of 16GB in an FAT-32 file system only. File formats supported by Braille Me are ".TXT" and ".BRF" encoded in 'ANSI (ASCII)'.

- **Does it support text to braille conversion?**

There is no on-board software to cater to this need at the moment. However, we are planning to facilitate the translation in future versions of the device. As for now, we

have provided the 'Braille Me translator', a desktop application to convert your files into any of the compatible formats and then transfer them onto the SD card for use.

- **How is Braille Me so low cost as compared to other braille displays?**

Most other braille displays use piezo-electric technology for their braille cells, which is an expensive approach. Innovision has developed and patented a unique magnetic framework which reduces cost immensely while keeping the functionality intact. Also, cost optimization has been implemented at every stage from manufacturing to distribution to cater to a larger customer segment.

- **Does it have audio input or output text to speech?**

No, audio and speech recognition hasn't been incorporated into the device. Braille Me is specifically designed to cater to a person's braille requirements.

- **With which screen readers does the device work?**

We have provided support for NVDA for Windows, BrailleBack for Android and Voiceover (currently in development) for MacOS and iOS.

- **How many bluetooth devices can be paired with Braille Me?**

Braille Me supports a maximum of four paired devices at a time.

2. First time user

- **What is the proper orientation of the device?**

Place the device on a flat surface in front of you with the row of braille cells farther away from you and the sloping edge nearer to you. This is the proper operational orientation.

- **Why is the device 6 dot and not 8 dot? How do you show cursor in 6 dot?**

The pressing dilemma at hand was whether to go for an 8 dot system for its utility or a 6 dot device for its cost efficacy. Therefore, we came up with a solution to integrate the cursor in the 6 dot module, thus, keeping the device operational and efficient at a lower cost while also being more accessible to the public.

There are two modes of displaying the cursor: In the first mode, the cursor blinks on the display, with the pins of the character being highlighted upright while the rest of the pins in periodic up and down motion. In the second mode, the manually triggered mode that is, the user has to press the Up + Down navigation keys simultaneously which results in a blank screen with only the cell accommodating the cursor to show up. Press the same combination keys again to revert changes. Refer to the manual for more details.

- **How long does the product work on batteries before it needs a recharge?**

The product takes about 1.5 hours to charge completely and following that, it can be used continuously for 10-12 hours before another recharge.

- **How do I pair a device via Bluetooth?**

First go to 'Bluetooth; option in home screen, then go to 'pair a new device' option. As you enter this menu item Braille Me will show a message – 'waiting for pair request'. Initiate pairing from the smartphone, and you will hear a sound once pairing is completed on the smartphone. Your smartphone might ask to confirm the pairing, please go ahead and press OK.

Wait for a few secs until Braille Me changes its message from "waiting for pair request" to an edit box asking "Enter paired device name." Enter a short name which you can use to identify your smartphone.

Once you receive a notification on Braille Me saying – "Paired to device name," the pairing process is complete.

- **Where can I find the product demonstration videos?**

To view the demo videos, follow this URL link:

<https://www.youtube.com/channel/UCmkFbftI7MKUS2Pjfx4ppg/videos>

- **Where can I find Braille books online?**

1. Sugamya Pustakalaya

<https://library.daisyindia.org/NALP/welcomeLink.action>

2. NLS

<https://www.loc.gov/nls/>

3. Bookshare Braille Central

<https://www.bookshare.org/cms/help-center/training-and-resources/braille>

- **Can I use my Kindle app on my Android smartphone?**

Make sure your BrailleBack or Voiceover is running in the background. Now on the Amazon Kindle app, browse for the book you are looking for. You can either pay for a subscription or read any of the free books available. Once a certain page is open, the screen reader reads the contents of the page and reproduces it on the Braille Me display.

- **How do I clean the braille display if needed?**

Brushing away the display is a good option. Might also use a soft piece of cloth and press gently. Strictly avoid using water/Colin and other related wet media for cleaning purposes.

- **How do I read a document in my native language on Braille Me?**

You can simply run the Braille Me translator software which converts the file into any of the supported formats. Once converted, copy it into your SD card and plug it into Braille Me. Now you can browse to the file and start reading it.

3. Existing User

- **What are the warranty plans for the product and/or any of its accessories?**

We provide a 1 year limited warranty for all accessories including the Braille Me unit itself.

- **My Braille Me has trouble pairing with my mobile device?**

One of the various possible issues may be that there are already 4 active devices paired to your Braille Me or you are trying to pair a device with a name that had already been reserved for some other device. To rectify the situation, you can delete a paired device from your Braille Me to accommodate the new device and try a different name for the former and latter situations respectively. For other issues refer to the user manual or contact customer support.

- **My computer doesn't recognize Braille Me as mass media storage?**

This may be due to the fact that the user hasn't followed the proper procedure to connect the device with the desktop. As long as 'connected as a mass device' is not shown on the display, the device will not connect to the desktop. If even after the message is displayed, there is trouble over connection, try unplugging the USB cable and re-plugging it. If the problem persists, contact customer support.

- **One of my typing keys has stuck / There's a problem with one or more pins in the display panel?**

If a key is jammed, an error message is displayed on the screen indicating the user about the issue. You can try prying free the key but do not try too hard. If any of the dots are stuck try cleaning it using a soft brush or cloth. If the issue persists then contact customer support.